

8. Sinogram/Fistulagram Pelvic Protocol

Sino gram/Fistula gram Pelvic Protocol

This need to be supervise by facilitator/Second Call

3 PL localizer free breathing

T2 TSE fs coronal

T2 TSE fs Axial

T2 TSE fs sagittal

T2 SPACE axial iso (**Call Radiologist to position SPACE and include entire anus and perianal tissues in FOV**)

T1 Vibe Dixon fs axial (4 points)

Diffusion b50, 800 BH (Match slice positioning with SPACE)

T1 VIBE fs axial breath hold

T1 VIBE fs sagittal breath hold

T1 VIBE fs coronal breath hold

Contrast (include entire anus and perianal tissues in FOV)

T1 VIBE fs axial +c (50s)

T1 VIBE fs axial +c (90s)

T1 VIBE fs sagittal+c (2 min)

T1 VIBE fs coronal+c (3 min)

T1 TSE axial +c (4 min)