

## 1. Routine Liver with MRCP with Contrast

<b>Routine Liver with MRCP with Contrast</b>
3 Plane Localizer free breathing
3 Plane Localizer breath hold
T2 Haste Coronal
T2 Haste axial
T1 ViBE Dixon breath hold /Fat quantification
T1 VIBE fs axial breath hold
<b>Post contrast</b>
T1 VIBE fs axial 30s
T1 VIBE fs axial 70s
T2 Haste axial fs
T1 VIBE fs axial 3 min
T2 TrueFISP axial bh
T1 VIBE fs axial 5 min
T1 Vibe coronal
Diffusion 50, 800
T2 Haste fs coronal thick slab breath hold