

Neck MRI

- Coronal T1
- Axial T1
- Axial T2 FS
- Coronal STIR
- Sagittal T1
- Axial DWI/ADC

- Post-contrast:
 - Coronal, axial and sagittal T1 FS

- DO NOT FAT SUPPRESS THE PRE-CONTRAST T1 SEQUENCES! These should be weighted so that normal fatty marrow is hyperintense. On the Skyra, please ensure that the T1-weighted sequences are performed as TSE/FSE T1.
- Can discuss with neuroradiologist to select coverage region